

**GO-GO<sup>®</sup> ES 8**

**OUTSTANDING VALUE.  
EXCEPTIONAL DURABILITY.**

## GO-GO<sup>®</sup> ES 8

The **Go-Go<sup>®</sup> ES 8** Scooter from Pride<sup>®</sup> delivers high-performance at an affordable price. The lightweight Go-Go ES 8 features easy disassembly, which permits simple frame separation with only one hand. Designed for exceptional durability, the Go-Go ES 8 provides reliability at an outstanding cost.



**Pride<sup>®</sup>**  
Mobility Products Ltd.  
— LIVE YOUR BEST<sup>®</sup> —

(Tel) 01869 324 600 • (Fax) 01869 323 070 • [www.pride-mobility.co.uk](http://www.pride-mobility.co.uk)

# Go-Go® ES 8

Weight capacity:	19 Stone 9 lbs. (125 kg)
Maximum speed <sup>1</sup> :	Up to 4 mph (6.4 km/h)
Ground clearance <sup>2</sup> :	5 cm (2") at mid-frame
Turning radius:	108.6 cm (42.75")
Overall length <sup>2</sup> :	102.2 cm (40.25")
Overall width <sup>2,6</sup> :	49.53 cm (19.5")
Front tyres:	5 cm x 17.8 cm (2" x 7") solid
Rear tyres:	6.35 cm x 20.32 cm (2.5" x 8") solid
Total weight without batteries <sup>2,5</sup> :	35 kg (77 lbs.)
Heaviest piece when disassembled <sup>2</sup> :	14.5 kg (32 lbs.) front section
Standard seating:	Type: Foldable, compact, vinyl-covered molded plastic Weight: 10 kg (22 lbs.) Material: Black Vinyl Dimensions: 45.7 cm (18") overall width 43.2 cm (17") effective depth
Drive system:	24-volt DC motor, sealed transaxle, rear-wheel drive
Dual braking system:	Regenerative and electromechanical
Battery requirements <sup>3</sup> :	Two 12-volt, deep cycle Size: 12 AH (standard) Weight <sup>4</sup> : 12 AH: 7.3 kg (16 lbs.) battery pack
Battery charger:	1.5 amp, off-board

Colours:



Red



Blue

## Features

- Front frame-mounted seat post offers maximum stability
- Easily disassembles into 6 lightweight pieces for convenient transport and storage
- Modular design for easy serviceability



Easy disassembly

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. This specification can be subject to a variance of +10%, -5%.

<sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

<sup>3</sup> AGM or gel cell type required.

<sup>4</sup> Battery weight may vary based on manufacturer.

<sup>5</sup> Includes standard seat.

<sup>6</sup> Overall scooter width at armrests is 59.7 cm (23.5").

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL  
(Tel) 01869 324 600 • (Fax) 01869 323 070 • [www.pridemobility.co.uk](http://www.pridemobility.co.uk)



DEALER INFO HERE

**Pride**<sup>®</sup>  
Mobility Products Ltd.

UK\_SC\_GO-GO ES8\_SS\_2-1-13



**PRIDE**<sup>®</sup>  
PERFORMANCE MOBILITY  
— LIVE YOUR BEST<sup>®</sup> —