

The new B400


Manoeuvrable and compact – for indoor and outdoor use

Quality for life



What did the B400 power wheelchair stand for in the past? A compact and stable construction: thanks to this power wheelchair the user could easily pass through narrow doorways and manoeuvre under tables and workstations.

What does the B400 stand for now? The same things. But also some new ones: the B400 has a new look and a modified suspension for improved comfort.



Mobile in all situations

A wheelchair should be small and ideal for all daily activities – whether at home, in a restaurant or on the go.

Compact construction

With modified design and added comfort

**THE NEW B400 WILL BE AVAILABLE
FROM SPRING 2013**

The B400 features a compact design with an overall width of 58 cm. Thus it can easily pass through narrow doorways. The low seat height of 40 cm allows the user to manoeuvre the wheelchair under tables and workstations with no problems. The compact construction is the greatest advantage of the B400. It is also a good choice for users on the go. It can easily negotiate curbs and inclines.

The new standard suspension on the drive unit provides additional comfort on long journeys. The design has been changed to give the B400 a fresh new look.

The small, versatile B400 is ideal for everyday requirements, whether at home, whilst shopping or in a restaurant. It features a long-proven technology and outstanding quality. Compact, manoeuvrable and robust – the B400 makes sure the user is always on the right track.

Improved functioning, improved appearance Innovations to the B400



**Frame shape:
optimised**



**Enhanced driving
comfort:
drive unit with
suspension support
in combination with
50Ah AGM batteries**

Now with the option '10 km/h'



Chassis:
now also available in
the colour 'cream'



Splash guard:
design adapted to
new frame shape



Fork:
new design



The new features

The B400 has an additional drive unit suspension which in combination with 50Ah AGM batteries is a standard feature. This allows users to travel long distances even on rough terrain in comfort. Another advantage for outdoor use: the B400 comes with an optional 10 km/h feature.

The new design

Cream with brown lettering or black with silver lettering? Your customer can select the appearance of his or her wheelchair. The frame was changed to give the B400 fine contours. The fork and the splash guard shapes are also more aesthetic.

Options

- Mechanically elevating footrests
- Electrically elevating footrests
- Curb climbing assist
- Lap belt
- Electric lighting
- Power back angle adjustment, 30°
- Caster swivel lock
- Luggage carrier
- Rear-view mirror
- Puncture-proof tyres
- Caster wheel suspension
- Attendant control
- Contoured seat
- 63 Ah (C5) gel batteries
- Control panel holder, swing-away in parallel
- Control panel holder, height-adjustable and swing-away in parallel
- Crutch holder
- Headrest

Colours



black



cream

Technical data

Top speed	6 km/h, 7.2 km/h or 10 km/h
Battery capacity	50Ah (C5) AGM batteries/63Ah (C5) gel batteries
Range	25/35 km
Battery charge time	max. 10 h
User weight	140 kg
Net weight	95 kg
Max. obstacle negotiation height	5 cm/10 cm (with curb climbing assist)
Climbing ability	12 %
Width	58 cm
Length including footplates	106 cm
Turning radius	80 cm
Seat height	40–57 cm
Seat width	38–48 cm
Seat depth	36–50 cm
Seat inclination	0° / 3° / 6°
Armrest height	24–36 cm
Arm pad length	26 cm
Back height	45/55 cm
Back width	40/44 cm
Back angles	-9° / 1° / 11° / 21° or 0° / 10° / 20° / 30°
Lower leg length	25–49 cm (standard seat), 28–54 cm (contoured seat)
Frame colour	cream or black
Control unit	VR2

Otto Bock HealthCare GmbH
Max-Näder-Straße 15, 37115 Duderstadt/Germany
T +49 5527 848 1590 · F +49 5527 848 1676
reha-export@ottobock.de · www.ottobock.com