

RASCAL
MOBILITY



Ultralite
480

Technical Specifications

Ultralite 480

Overall length	1020mm (40")
Overall width	510mm (20")
Overall height	850mm (33")
Maximum rider weight	18 stone (113kg)
Batteries	2 x 12 amp/hr
Ground clearance	5cm
Total weight	41kg (90lb)
Heaviest component weight	13kg (29lb)
Maximum range	Up to 8 miles
Maximum speed	4mph
Colour	Red / Black / White

Due to a policy of continuous improvement, Electric Mobility reserves the right to change specifications without prior notice.



0800 252 164 sales@electricmobility.co.uk electricmobility.co.uk



Please contact your local
Rascal Mobility Dealer

