

Invacare® | Pronto® M61™



A fresh approach to mobility

The sleek and stylish *Invacare Pronto M61* is an affordable centre wheel drive power chair with 6 wheels on the ground for optimum safety and stability. It features *Invacare's Sure Step®* technology for smoother driving over lower thresholds. Designed around a compact base, the *Pronto M61* provides excellent manoeuvrability especially in confined spaces, and the electric seat riser provides 120 mm of elevation.

INVACARE
Yes, you can.®



Features and Options

Optimum Comfort

The padded, shaped van seat ensures you are comfortable for day long use, while the popular configurable seating system offers a wide range of seat sizes allowing you to achieve an optimum fit. For easier transfers into or out of the **Pronto M61**, the seat will rotate 90 degrees in either direction, and the van seat armrest can flip up or alternatively the armrests can be removed altogether.



Sure Step® technology gives a smooth ride over low thresholds of up to 5cm.

Upwardly mobile

With its centre wheel drive and its compact dimensions, the **Pronto M61** gives excellent manoeuvring in very small spaces and has a very tight turning circle. The simple to operate seat will electrically raise you up to 120 mm. In the elevated position it provides increased accessibility, which is ideal in an office environment or when you are out and about.



A range of seating options are available including the **Laguna** seat and back style.



A range of angle adjustable legrests and footplates available. (Configured version only)

Technical data

Pronto M61 Standard	455 - 500 mm	410 - 460 mm	0 - 25°	890 - 940 mm
Pronto M61 Configured	330 - 430 mm	300 - 480 mm	12° (fixed)	960 - 1100 mm
Pronto M61	470 - 600 mm ¹ 500 - 630 mm ²	610 mm (excluding joystick)	680 mm ³ 990 - 1030 mm ⁴	R = 495 mm (front and rear turning radius)
Pronto M61	9° / 16 %	2 x 30 Ah	2 x 180 W	16 km 6.4 km/h
		Frame colours		
Pronto M61	70 kg (94 kg with batteries)	136 kg		
			Pearl Blue	Red Pearl

Invacare Limited
Pencoed Technology Park - Pencoed - Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail: ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2013 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare reserves the right to alter product specifications without prior consultation.

Pronto - UK - 04/2013 - 1468372



Yes, you can.®

- 1/ Seat height is measured including cushion, with 4° tilt angle
- 2/ Seat height is measured including cushion, with 9° tilt angle
- 3/ with footboard up
- 4/ with footboard down, with side mounted legrests