

Neon

First class travel...



Neon

The Mini-Max Principle

Compact, folding and lightweight for easy stowage in a plane, train or car.

Lightweight frame weighs just 7.1kgs with intuitive folding cross brace and retaining strap, removable legrests and fold down backrest.



Neon Fixed Front version



Neon Swing-Away version

Getting free and independent



First class travel through fine adjustment... journey after journey.

Infinite fine tuning possibilities to ensure perfect driving ergonomics can be achieved and maintained.

Fashionable style and looks in a combination to suit you.

Up to 3 million option combinations with over 20 design choices.

3.7 million colour combinations including six upholstery colours and 27 brilliant and matt paint finishes for your frame, hubs, forks, wheels, handrims and footplates.



Modern open frame design (Fixed Front)



Auto-folding footboard (option)



Small folding frame (Fixed Front)



Internal cross-brace ensures optimum rigidity and drive performance



Simple seat height adjustment with preset factory "centre of gravity" settings



Improved, durable folding backrest system offering angle adjustment from -12° to $+8^{\circ}$



Neon Swing-Away: with swing-in and swing-out function





Fast folding frame and back with quick release legrest ensures convenient transportation





Wide range of side guards and armrest options



There is room for the **Neon**
even in the smallest car.



Swing-away:



Fixed Front:



Technical data sheet

27 individual frame colours

Folded width: 37cm (14.5") (SA) / 38cm (15") (FF)

6 individual sling colours

Optional seat depth growth 2.5cm (1") and 5cm (2")

Max. user weight: 115kg* (18 stone) / 140kg** (22 stone)

Weight of the heaviest part: 7.4kg (16.3lbs) (SA) / 7.5kg (16.5lbs) (FF)

Total weight: From 9.9kg (21.8lbs)

Transit tested to ISO 7176 Part 19 and passed: Yes***

Find out more today!

Visit www.SunriseMedical.com for more information, to find your local dealer or to request a demonstration.

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.



Get SOCIAL with QUICKIE

Keep up-to-date with the latest news on the Quickie Neon and other Quickie products through our Facebook page:

www.facebook.com/QuickieWheelchairsUK



Sunrise Medical Limited
Thorns Road
Brierley Hill
West Midlands
DY5 2BQ

Tel: +44 (0) 845 605 66 88

Fax: +44 (0) 845 605 66 89

www.SunriseMedical.com

* Seat width 42cm (16.5") ** Above seat width 42cm (16.5") *** Quickie Neon Swing-Away model only.