

NEW



**LiteLock™**  
system

(patent pending)



## The kind of Scooter you want to drive

The new *Invacare Colibri* is a stylish and colourful boot scooter, designed for those who enjoy their independence. Its flexibility means the scooter fits around your lifestyle whilst remaining effortless to use. The *Invacare Colibri* is all about choices. So depending on where you go, how far you travel and what your favourite colour is, **Colibri** has something for you.

Thanks to the unique *Invacare LiteLock™* system the scooter can easily be taken apart without the need for tools and fits neatly into any car boot.

You can choose from a range of modern colour panels to match your style.



**Yes, you can.®**



**Features and Options**

**LiteLock system - easy to take apart and reassemble**

The ingenious and unique *Invacare LiteLock™* system automatically releases once the battery box has been removed. It takes less than a minute to take the scooter apart – without tools. Its front section can be simply lifted and placed in the boot of a car, along with the other sections. It's just as easy to reassemble the scooter when you need it again.



**Effortless dismantling**

Don't think - do it! Less than a minute is all that's needed to disassemble the **Colibri** for transportation.



**Option of a bigger wheel set and battery**

Perfect for travelling longer distances, on bumpier terrain, such as cobbled streets. Choose where to charge your battery: directly on the **Colibri** or detached from it.



**Increased comfort**

The soft, padded seat provides comfort during driving. For increased comfort there is also the optional 'easy fit' seat suspension kit.

**Technical data**

<i>Invacare Colibri</i> 465 mm	400 mm	370 - 420 mm	385 mm	505 mm	1010 mm
<i>Invacare Colibri</i> 44 kg (12 Ah batteries) 49 kg (18 Ah batteries)	136 kg	Weight of heaviest part 15.6 kg	2 x 12 Ah 2 X 18 Ah	2200 mm	
<i>Invacare Colibri</i> 45 mm	6° / 10.5%	200 W	Up to 16 km	6 km/h	

Invacare Limited  
Pencoed Technology Park - Pencoed -  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS - County  
Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

© 2014 Invacare International Sàrl  
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Colibri - UK - 03/2014 - 1576204


**Shroud colours**



Sapphire Blue

Jasper Red

Jade Green



Sunstone Orange

Marble White



**Yes, you can.®**