

RASCAL
MOBILITY



Pioneer

Technical Specifications

	Pioneer	Frontier
Overall length	1470mm (57.9")	1430mm (56.3")
Overall width	645mm (25")	650mm (26")
Overall height	1300mm (51")	1020mm (40")
Maximum rider weight	25 stone (159kg)	21 stone (136kg)
Batteries	2 x 75 amp/hr	2 x 55 amp/hr
Ground clearance	9cm	9cm
Total weight	130kg	99.5kg
Maximum range	Up to 35 miles	Up to 25 miles
Maximum speed	8mph	8mph
Wheel size (front)	12"	10"
Wheel size (rear)	12"	11"
Turning radius	1800mm	1400mm
Colour	Red / Black	Red / Black

Due to a policy of continuous improvement, Electric Mobility reserves the right to change specifications without prior notice.



0800 252 164 sales@electricmobility.co.uk electricmobility.co.uk



Please contact your local
Rascal Mobility Dealer

