



Yes, you can.®

➤ **Invacare® Action³_{NG} Comfort**

For additional comfort and positioning

The *Invacare® Action³_{NG} Comfort* now offers optimal lateral positioning and additional seating comfort with its new smooth Dartex® fabric which is hygienic and easy to maintain. It is the new comfort folding wheelchair.



➤ New comfort

The **Invacare® Action^{®3} NG Comfort** offers a 30° reclining backrest fitted as standard, with seat and backrest anatomic preformed soft cushions. In addition, the **Action^{®3} NG Comfort** now features height and depth adjustable armrests, and legrests with knee protection to provide elderly people with optimal comfort.



➤ Dartex® ergonomic seating

The ergonomic seating (seat and backrest cushions, armrests and headrest can be quickly removed) is covered in Dartex®, a supple and very smooth coating, thus offering unmatched positioning and comfort.

Standard features



Legrests

The multi-adjustable legrests fitted with articulated footplates provide perfect positioning.



Headrest


Multi adjustable; comfort and positioning are optimal



Armrests

The combinations of adjustments ensure ultimate seating positioning.

Technical data




40,5 / 43 / 45,5 / 50,5 cm




45 cm



53* / 50* / 48* cm



56 cm




17 - 26 cm




37 - 49 cm



0° - 30°



59,5 / 62 / 64,5 / 69,5 cm




108 cm




113 cm



86 cm



32 cm



26,5 kg



125 kg



Crashtested ISO 7176-19**


Frame colour



Slate Grey



Ice Brown



Ocean Blue



Anthracite Black



Dark Red

Upholstery



Grey Dartex

* Seat height including ergonomic cushion (6 cm thickness)

** Our wheelchairs comply above ISO standards and with ISO 10542, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.



INVACARE® Limited
Pencoed Technology Park · Pencoed · Bridgend CF35 5HZ
Tel: +44 | 656 776 222 · Fax: +44 | 656 776 220
E-mail: uk@invacare.com · www.invacare.co.uk
Sales Order E-mail: ordersuk@invacare.com

INVACARE® Ireland Ltd
Unit 5 · Seatown Business Campus
Seatown Road · SWORDS · County Dublin
Tel. +353 | 8107084 · Fax +353 | 8107085
Email: Ireland@invacare.com · www.invacare.ie

A3 Comfort · UK · 01/10 · Part-No. 1504618
© Invacare® International Sàrl 2004-2010. All rights reserved.

