

Jazzy[®]

SELECT 6

The Jazzy Select[®] 6 provides exceptional front stability combined with the advantages of in-line motor technology. Active-Trac[®] Suspension and 6" front and rear casters offer excellent performance indoors and out. An impressive selection of standard convenience features make it the easiest power chair to use.

Jazzy Select[®] 6
shown with basic Jazzy[®]
Armor[®] (standard)



Pride[®]
Mobility Products Corp.

Specifications

Weight capacity	300 lbs.
Maximum speed	Up to 4.2 mph
Ground clearance	1.75"
Turning radius	21.5"
Overall length	36.625"
Overall width	23.5"
Drive wheels	10" solid
Front wheels	6" solid
Rear wheels	6" solid
Suspension	Active-Trac®
Drivetrain	Two-motor, in-line, mid-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	40 amp, PG GC 3 controller
Specialty controls	N/A
Battery charger	Off-board, 3.5 amp
Per-charge range	Up to 11 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	U-1
Battery weight	25 lbs. each
Base weight	100 lbs.
Standard seat weight	37 lbs. (high-back)
Max. seat dimensions	20" x 20"
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

Models

Jazzy Select 6-SS (Code: K0822)
Jazzy Select 6-C (Code: K0823)

Features

- 6" casters on front & rear for maximum stability
- In-line motor technology for increased efficiency, torque and range
- Articulating rear caster beam
- Easy access charger port and circuit breaker
- Height and depth-adjustable foot platform

Options

- Power elevating seat
- Depth-adjustable, fold-flat seat
- High-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

Recommended Lifts

- Backpacker® MV
- Backpacker® Plus
- Outlander
- PrideBoom 250



Color Choices



Red



Blue

Pride®
Mobility Products Corp.
LIVE YOUR BEST®

(US) 800-800-8586 • (Canada) 888-570-1113 • www.pridemobility.com