

SHARK RS

QUICKIE
Live without limits.



Designed for champions.
Made by Quickie.

>>>> Introducing the new Quickie Shark RS - the choice of champions

The Shark RS's new innovative frame is designed to allow for the best individual positioning. Enclosing the driver, when force is introduced, it's converted directly into propulsion. The low centre of gravity and wide wheelbase enable higher cornering speeds whilst maintaining maximum stability.

- **NEW** innovative frame design
- **NEW** lower weight
- **NEW** design improvements
- **INCREASED** free-moving space

>>>> A proven winning design



On Sunday 9th October 2011, Vico Merklein (pictured left) showcased the stealthy credentials of the new Quickie Shark RS handbike by taking victory in the 11th Handbike marathon in Heidelberg, Germany.



Practice>Determination>The NEW Quickie RS handbike>A winning combination>>>>

>>>> Performance and style

The NEW Quickie Shark RS sets a new standard in driving performance.

Weight reduction combined with design improvements.
Best design in class.



An innovative and ergonomic backrest combines independent and adjustable elements for enhanced performance.



Uses NEW DT-SWISS carbon wheels to give unique performance and excellent design.






A NEW rear axle with forged parts combines high stability with low weight.

The NEW front fork offers a lightweight yet rigid solution.

>>>> Technical Specifications

SHARK RS



- WEIGHT:** from 12.9kg
- FRAME WIDTH:** 34 - 42cm (13.5 - 16.5") in 2cm (0.8") increments
- REAR:** angle adjustable 20° to 60°
- CRANK LENGTH:** 15.5 - 17.5cm (6 - 7") in 0.5cm (0.2") increments
- GEAR:** 3x10 Shimano
- BRAKE:** SRAM
- COLOURS:** 30 colours available

Customisation available through



Ask your dealer for more details.



Sunrise Medical Limited

High Street
Wollaston
Stourbridge
West Midlands
DY8 4PS
Tel.: +44 (0) 845 605 6688
Fax: +44 (0) 845 605 6689
www.QuickieWheelchairs.eu/Shark
www.facebook.com/QuickieWheelchairsUK

Find out more today!

Visit www.QuickieWheelchairs.eu/Shark
for more information, to find your local dealer
or to request a demonstration